


# TEOAE

## Transient-Evoked Otoacoustic Emission Module



Normal-hearing adult showing a Pass result in 5 bands. Green patch shows the TEOAE spectrum. Orange patch indicates the noise floor.

### TE: Transient-Evoked Otoacoustic Emission

Non-invasively and objectively measure inner-ear function.

- Safe and effective for newborns to adults
- Takes just seconds
- Allows you to utilize factory preset protocols or design your own
- For use in hospitals, clinics, private practice, and research institutions
- Instrumentation employs:
  - A series of stimuli (clicks or chirps)
  - In-the-ear probe calibration
  - Greater control of test parameters in research mode

Combine with MEPA3 wideband immittance software for a complete middle- and inner-ear diagnostic system.



# HearID 5 + TE Module - Technical Specifications

Transient-Evoked Otoacoustic Emission System

## HearID + TE Specifications

To run HearID + TE, you need a HearID software and hardware system with the TE Module.

- Clinical and Research modes.
- Preset measurement protocols included.
- Synchronized spontaneous otoacoustic emissions mode.
- Stimulus Types: rectangular click, biphasic rectangular click, Gaussian click, differentiated Gaussian click, flat group-delay chirp, linear group-delay chirp, Dau group-delay chirp, Shera group-delay chirp.
- Customizable stimulus bandwidth, duration, level, white or pink spectrum, polarity.
- Adjustable inter-stimulus interval.
- Linear and nonlinear stimulus presentation.
- Adjustable number of passes and maximum rejected passes.
- Adjustable TEOAE response filters.
- Adjustable response window.
- Adjustable graphical noise rejection level, now with option to automate.
- Observe stimuli acoustically and electrically in real time.
- In-the-ear probe calibration, including option for spectrum compensation.
- Pass / Refer results clearly indicated.
- Protocols can be customized and saved.
- Improved stopping rules and pass/refer criteria settings.

## Environmental Requirements

For valid measurements:

- The air temperature must be within +59°F to +95°F (+15°C to +35°C).
- Relative air humidity must be between 30% and 90%.
- Ambient noise should be minimal; however, a sound-treated booth or room is not necessary.

## Regulatory

HearID TE is a legally marketed device (FDA 510(k)) and is CE marked.



## HearID + TE Hardware

- 24-bit proprietary USB Audio Processing Unit.
- ER10C sound probe from Etymotic Research.
- Disposable foam and rubber eartips for comfortable, stable probe insertion. Seven sizes available, suitable for infants to adults.
- Cutoff 2-cc syringe, for use as a test cavity.
- Laptop computer not included but must meet the minimum requirements below.

## Computer System Minimum Requirements

- IBM™ compatible laptop with a Pentium™ class processor at least 1.80GHz.
- Windows XP™ (SP3, 32-bit), Windows Vista™ (32- or 64-bit), or Windows 7™ (32- or 64-bit) operating system.
- For XP systems: at least 512MB RAM and at least 20GB hard drive.
- For Vista systems: at least 1GB RAM and at least 50GB hard drive.
- For Windows 7 systems: at least 2GB RAM and at least 80GB hard drive.
- One available USB 2.0 port.
- Recent version of Adobe Reader (for displaying the electronic manuals).
- EN60950 (IEC950)-compliant power supply.

## Order Information

Product Identifier: HearID+TE

## Contact Mimosa Acoustics

### Mimosa Acoustics, Inc.

335 N. Fremont St.  
Champaign, IL 61820, USA  
Phone +1-217-359-9740  
Fax: +1-217-359-9741  
[sales@mimosaacoustics.com](mailto:sales@mimosaacoustics.com)  
[www.mimosaacoustics.com](http://www.mimosaacoustics.com)

 [facebook.com/MimosaAcoustics](https://facebook.com/MimosaAcoustics)

 [@MimosaAcoustics](https://twitter.com/MimosaAcoustics)

